

**SemiVolatile Organic Compounds in DW by Liquid/Solid Extraction and GC/MS
EPA 525.2, Revision 2.0****Page 1 of 3**

Facility Name: _____ EPA ID _____

Assessor Name: _____ Analyst Name: _____ Inspection Date _____

Method Specific Requirements**Reference****Y****N****N/A****Comments**

Records Examined: _____

Date of Analysis: _____ Date(s) of Sample Preparation: _____ Analyst: _____

Is a series of 6 calibration standards prepared for all compounds of interest? (If all method analytes are to be determined, 2 or 3 sets of calibration standards will likely be required). *[response factor (if used) <30%RSD]

5.10.2.a
**EPA 7.10,
10.2***

Do at least 5 calibration points bracket the expected analyte concentration range?

5.10.2.a
EPA10.2.5.1

Are samples extracted within 14 days of collection and the extracts stored at 4°C and analyzed within 30 days of collection (except for carboxin, diazinon, disulfoton, disulfoton sulfoxide, fenamiphos, and terbufos, which are extracted immediately)?

5.10.2.a
EPA 8.3

Is the recovery of the Internal Standard in laboratory blanks and samples in excess of 70%? (Absolute peak areas for internal standard and surrogates must also be checked)

5.10.2.a
EPA 9.4

Is the laboratory fortified blank (LFB) containing each analyte analyzed with each batch of samples as a group within a 12 hour work shift?

5.10.2.a
EPA 9.6

Is a continuing calibration check analyzed each day or at the beginning of each period in which analyses are performed not to exceed 12 hours.

EPA 10.1

Is the tune of the GC/MS checked each day with DFTPP?

5.10.2.a
EPA10.2.2

Is the breakdown of Endrin and DDT checked each day and maintenance performed if the breakdown of either exceeds 20%? (GC injector port Degradation <20%)

5.10.2.a
EPA10.2.2

Is the medium calibration standard check to ensure that anthracene and phenanthrene are separated by baseline? Benz(a)anthracene and chrysene are separated by a valley whose height is less than 25% average peak height of these two? (every 12 Hours)

5.10.2.a
EPA10.2.4.1NOTES/COMMENTS:

SemiVolatile Organic Compounds in DW by Liquid/Solid Extraction and GC/MS EPA 525.2, Revision 2.0					Page 2 of 3
Method Specific Requirements	Reference	Y	N	N/A	Comments
Is the system able to recognize at least 99% of the compounds in the medium calibration standard?	5.10.2.a EPA10.2.4.2				
If the RSD of any analyte or surrogate mean RF exceeds 30%, are linear regression curves used?	5.10.2.a EPA10.2.6.1				
Is an external source quality control sample analyzed at least quarterly?	EPA 9.9				
Is the Internal Standard and surrogate areas not allowed to decrease from the latest continuing calibration check by more than 30% or the initial calibration by more than 50%?	5.10.2.a EPA 10.3.4				
[] Is the Reference Factor (RF) for each analyte and surrogate within 30% of the mean value measured in the initial calibration? or [] Alternatively, if the linear regression is used, is the recovery amount for each analyte +/- 30% of the true value? (A single point calibration can be used if there are hits for a few analytes that did not meet criteria.)	5.10.2.a EPA10.3.5				
Sample Preparation:					
Analyst: _____		Date: _____			
Do records indicate that the sample has been checked for residual chlorine and that the pH is less than 2?	5.12.3.1.a EPA 8.2				
Cartridge Extraction:					
After passing the sample through the LSE cartridge, is the separatory funnel and sample container rinsed with ethyl acetate and methylene chloride used to elute the cartridge?	5.10.2.a EPA 11.1.2				
Is the extract concentrated in a warm water bath under a gentle stream of nitrogen to a volume greater than 0.5 mL?	5.10.2.a EPA 11.1.5				
Disk Extraction:					
After passing the sample through the disk, is the sample container rinsed with ethyl acetate which is allowed to elute through the disk?	5.10.2.a EPA 11.2.6				
Are two aliquots of 1:1 ethyl acetate: methylene chloride allowed to elute through the disk?	5.10.2.a EPA 11.2.8				
<u>NOTES/COMMENTS:</u> 					

**SemiVolatile Organic Compounds in DW by Liquid/Solid Extraction and GC/MS
EPA 525.2, Revision 2.0**

Page 3 of 3

Method Specific Requirements	Reference	Y	N	N/A	Comments
Is the extract concentrated in a warm water bath to between 0.5 mL and 1 mL?	5.10.2.a EPA 11.2.9				
Do the records include the identification of any standards, surrogates, internal standards, solvents, etc. used in the analysis and preparation of the samples?	5.12.3.1.f (Extra Info.)				
<u>NOTES/COMMENTS:</u>					

Extractable Organics Gas Chromatography Mass Spectrometry (GC/MS)

Method Specific Requirements

-Required Reagents & Standards-

Reverse-phase C-18 solid-phase disks or cartridges

Ethyl Acetate, Methylene Chloride, Methanol solid-phase conditioning reagents

Ethyl Acetate, the Methylene Chloride, as eluting solvents

Sodium Sulfate drying reagents

Endrin & DDT for GC degradation check

Internal Standards

DFTPP (Decafluorotriphenylphosphine) for tuning MS (Daily)