Alkalinity by Flow Injection Colorimetry- LACHAT QuickChem Method 10-303-31-1-A					
Facility Name:	VELAP ID				
Assessor Name:Analyst Name:	Inspection Date				
Relevant Aspect of Standards	Method Reference	Υ	N	N/A	Comments
Records Examined: SOP Number/ Revision/ Date	ds Examined: SOP Number/ Revision/ Date Analyst:				
Sample ID: Date of Sample Prepa	aration:		_ Da	ate of A	Analysis:
Are samples preserved and cooled to ≤6°C at the time of collection?	40 CFR Part 136.3				
Are samples analyzed within 14 days of collection?	40 CFR Part 136.3				
Was an initial MDL established?	9.2.4				
Is reagent water CO2 free (boiled for 5 min)?	7.1				
Is KHP Buffer prepared 5 g KHP/1L and pH adjusted to 3.1 +/- 0.05?	7.1				
Is a Zero concentration standard prepared and analyzed as part of the standard curve?	7.2				
Is a LRB analyzed with each batch and determined to be less than the MDL?	9.3.1				
Is a LCS analyzed with each batch and the recovery within +/- 10%?	9.3.3				
Is an IPC analyzed with each batch and the recovery within +/- 10%?	9.3.4				
Are a minimum of 10 percent of all samples (one per batch of 10) analyzed in duplicate and RPD results meet current lab acceptance criteria?	9.3.5				
Are a minimum of 10 percent of all samples (one per batch of 10) spiked and is the recovery within ±10%?	9.4.1				
Notes/Comments:					