

FLAME ATOMIC ABSORPTION METALS (Nitrous Oxide/Acetylene Flame)
SM 3111 D – 1999 (2011)

Facility Name: _____ VELAP ID _____

Assessor Name: _____ Analyst Name: _____ Inspection Date _____

Records Examined: SOP Number/ Revision/ Date _____ Analyst: _____

Sample ID: _____ Date of Sample Preparation: _____ Date of Analysis: _____

Relevant Aspect of Standards	Reference	Y	N	N/A	Comments
1) If samples require digestion, are standards and method blanks also subjected to digestion? (Recommend)	SM3111.A.5				
2) Does the calibration curve include at least three concentrations of standard metal solutions?	SM3111D.4.c				
3) Are standards prepared by appropriate dilution of stock metal solutions with water containing 1.5 mL conc. HNO ₃ /L?	SM3111D.3.k				
4) Is an additional standard (CCV) analyzed after every batch of 10 or fewer samples? (Recommended concentrations and criteria are in Table 3111:III.)	SM3111.A.7				
5) Is a matrix spike performed every 10 samples or each batch with an acceptable recovery of 85-115%?	SM3111.A.7				
6) Are hollow cathode lamps or electrode-less discharge (EDL) lamps used?	SM3111A.6.d				
7) Are fuel supplies are maintained at pressures slightly higher than controlled operating pressure of the instrument by suitable valves?	SM3111A.6.e				
8) Fuel & oxidant: [] Commercial grade acetylene, replaced when pressure has fallen to 689 kPa (100psi) acetylene; [] Acetylene cylinders are changed before pressure drops below 100 psi; [] NO copper or brass regulators, tubing, or fittings with >65% copper are used with acetylene; [] Nitrous oxide supply is equipped with a non-freezable regulator, or a heating coil is wrapped around an ordinary regulator to prevent flashback at the burner, or the instrument has automatic gas control system.	SM3111D.3.a SM3111D.3.b SM3111D.3.h.				

Notes/Comments

