

## Carbapenem Resistant Enterobacteriaceae (CRE) and Pseudomonas aeruginosa (CRPA) Testing Instructions

### I. CURRENT DCLS LABORATORY TESTING

Test Type	Method
Identification Confirmation	MALDI-TOF (Bruker Biotyper)
Antimicrobial Susceptibility Testing (AST)*	Broth Microdilution/Sensititre
Phenotypic testing - Carbapenemase Production	Modified Carbapenem Inactivation Method (mCIM)
Resistance Mechanism Identification	Real-Time PCR

**\* Submitting laboratories are REQUIRED to provide AST results when submitting suspected CRE/CRPA isolates to DCLS for testing. These AST results are assessed upon receipt to determine potential pan-resistance.**

#### TESTING METHODS

- **mCIM** – Phenotypic testing for carbapenemase production will be performed for all CRE and CRPA isolates.
- **AST** – The following drug Minimum Inhibitory Concentration (MIC) values will be determined for all CRE isolates and mCIM-positive CRPA isolates:
 

<ol style="list-style-type: none"> <li>1. Amikacin</li> <li>2. Aztreonam*</li> <li>3. Cefepime*</li> <li>4. Cefotaxime* (CRE only)</li> <li>5. Ceftazidime*</li> <li>6. Ciprofloxacin</li> <li>7. Colistin*</li> <li>8. Ertapenem (CRE only)</li> </ol>	<ol style="list-style-type: none"> <li>9. Imipenem*</li> <li>10. Levofloxacin</li> <li>11. Meropenem*</li> <li>12. Piperacillin/Tazobactam**</li> <li>13. Polymyxin B</li> <li>14. Tobramycin</li> <li>15. Trimethoprim/Sulfamethoxazole (CRE only)</li> </ol>
---	--

\* Testing required by the CDC's Antibiotic Resistance Laboratory Network (ARLN) for CRE and CRPA isolates.

\*\* Testing required by the CDC's ARLN for *Pseudomonas aeruginosa* isolates.

- **Real-time PCR** - Analyses will be performed on the following resistance marker genes for all CRE isolates and mCIM-positive CRPA isolates:
 

<ol style="list-style-type: none"> <li>1. KPC</li> <li>2. NDM</li> <li>3. OXA-48-like</li> </ol>	<ol style="list-style-type: none"> <li>4. VIM</li> <li>5. IMP</li> <li>6. <i>mcr-1</i> and <i>mcr-2</i></li> </ol>
--	--

### II. SPECIMEN TYPES

- Pure suspect CRE isolates on slant or plate media.
  - 1) Member of the Enterobacteriaceae family
  - 2) Resistance to at least 1 carbapenem:
    - imipenem, meropenem or doripenem = MIC of  $\geq 4 \mu\text{g/mL}$
    - ertapenem = MIC of  $\geq 2 \mu\text{g/mL}$

**Note:** *Morganella spp.*, *Proteus spp.*, and *Providencia spp.* have intrinsic elevated MICs to imipenem. MIC results for meropenem, doripenem, and/or ertapenem must be used to determine if these organisms should be forwarded to DCLS for testing.

- Pure suspect CRPA isolates on slant or plate media.
  - 1) Isolate identified as *Pseudomonas aeruginosa*
  - 2) Resistance to at least 1 carbapenem:
    - imipenem, meropenem or doripenem = MIC of  $\geq 8 \mu\text{g/mL}$
  - 3) Mucoid *Pseudomonas aeruginosa* isolates will NOT be tested

### **III. SPECIMEN TRANSPORT**

- Ship isolates at room temperature.
- Submit a completed DCLS Clinical Microbiology/Virology Test Request form with isolate, to include the following:
  - 1) Complete patient information (including address)
  - 2) Complete submitter information
  - 3) Date of specimen collection
  - 4) Specimen source
  - 5) Organism identification, denoting “suspected carbapenem resistance”
  - 6) Outbreak ID number (if applicable)
- Include copy of laboratory’s AST results in specimen package
- **Shipping Address:**  
Division of Consolidated Laboratory Services  
Attn: CRE Team  
600 North 5<sup>th</sup> Street  
Richmond, Virginia 23219

### **IV. RESULTS REPORTING**

- Turnaround time for reporting results will be within 6 business days of specimen receipt.
- Alert value results (**results requiring additional reference laboratory testing**) will be verbally reported within 1 working days of results. An alert email notification will be sent to the CDC as well.
- Non-alert-value, positive results (**results NOT requiring additional reference laboratory testing**) will be verbally reported to the submitting laboratory and VDH within 2 working days of results.
- Negative results will not be verbally reported.
- A hard copy report of final test results will be provided to the submitter by mail.
  - CRE final results will include: identification, AST, phenotypic, and molecular
  - CRPA final results will include: if mCIM negative, identification and phenotypic only; if mCIM positive, identification, AST, phenotypic, and molecular
- VDH will receive electronic reporting of results.
- Reference laboratory results for additional characterization testing will be provided to the submitting laboratory and VDH within 2 working days of receipt from the reference laboratory.

### **V. CONTACT INFORMATION**

- **For questions regarding testing, specimen collection or shipment, please contact one of the following:**
  - Dr. LaToya Griffin-Thomas (804-648-4480, ext. 281)
  - Microbial Reference Laboratory (804-648-4480, ext. 345 or 251)
  - Molecular Detection and Characterization Laboratory (804-648-4480, ext. 298)
  - DCLS 24/7 Emergency Mobile Number (804-335-4617)

### **VI. RESOURCES**

Updated information on CRE/CRPA testing can be found on the DCLS website, under the “Hot Topics” link:  
<https://dgs.virginia.gov/division-of-consolidated-laboratory-services/updates/hot-topics/>

Link to DCLS Clinical Microbiology/Virology Test Request Form:  
[https://dgs.virginia.gov/globalassets/business-units/dcls/documents/sample-kits/dcls-micro\\_viro-form\\_fillable\\_rw.pdf](https://dgs.virginia.gov/globalassets/business-units/dcls/documents/sample-kits/dcls-micro_viro-form_fillable_rw.pdf)