

## Carbapenem Resistant Organism (CRO) Testing Instructions

### • CURRENT DCLS LABORATORY TESTING

Test Type	Method
Identification Confirmation	MALDI-TOF (Bruker Biotyper)
Antimicrobial Susceptibility Testing (AST)*	Broth Microdilution/Sensititre
Phenotypic testing - Carbapenemase Production	Modified Carbapenem Inactivation Method (mCIM)
Resistance Mechanism Identification	Real-Time PCR

\* Submitting laboratories are **REQUIRED** to provide AST results when submitting suspected CRO isolates to DCLS for testing. These AST results are assessed upon receipt to determine suitability for testing and potential pan-resistance.

### TESTING METHODS

- **mCIM** – Phenotypic testing for carbapenemase production will be performed for all Carbapenem Resistant Enterobacteriales (CRE) and Carbapenem Resistant *Pseudomonas aeruginosa* (CRPA) isolates.
- **AST** – The following drug Minimum Inhibitory Concentration (MIC) values will be determined for all Carbapenem Resistant *Acinetobacter baumannii* (CRAB), as well as mCIM positive CRE and CRPA isolates:
 

1. Amikacin	10. Ciprofloxacin
2. Ampicillin (CRE only)	11. Ertapenem (CRE only)
3. Ampicillin/sulbactam (CRE/CRAB only)	12. Gentamicin (CRE only)
4. <b>Aztreonam* (CRE/CRPA only)</b>	13. <b>Imipenem*</b>
5. Cefazolin (CRE only)	14. Levofloxacin
6. <b>Cefepime*</b>	15. <b>Meropenem*</b>
7. <b>Ceftazidime*</b>	16. <b>Piperacillin/tazobactam**</b>
8. Ceftazidime/avibactam (CRE/CRPA only)	17. Tobramycin
9. Ceftolozane/tazobactam (CRE/CRPA only)	18. Tetracycline (CRE/CRAB only)
	19. Trimethoprim/sulfamethoxazole (CRE/CRAB only)

\* Testing required by the CDC's Antibiotic Resistance Laboratory Network (ARLN) for CRO isolates.

\*\* Testing required by the CDC's ARLN for *Pseudomonas aeruginosa* and *Acinetobacter baumannii* isolates.

- **Real-time PCR** - Analyses will be performed on the following resistance marker genes for CRAB isolates and mCIM positive CRE and CRPA isolates:

<u>CRE/CRPA</u>	<u>CRAB</u>
<ul style="list-style-type: none"> <li>• KPC</li> <li>• NDM</li> <li>• VIM</li> <li>• OXA-48-like</li> <li>• IMP</li> </ul>	<ul style="list-style-type: none"> <li>• KPC</li> <li>• NDM</li> <li>• VIM</li> <li>• OXA-48-like</li> <li>• IMP</li> <li>• OXA-24/40</li> <li>• OXA-23</li> <li>• OXA-58</li> </ul>

## • **SPECIMEN TYPES**

- Pure suspect CRE isolates on slant or plate media.
  - 1) Member of the Enterobacterales family
  - 2) Resistance to at least 1 carbapenem:
    - imipenem, meropenem or doripenem = **MIC of  $\geq 4$  µg/mL**
    - ertapenem = **MIC of  $\geq 2$  µg/mL**

**Note:** *Morganella spp.*, *Proteus spp.*, and *Providencia spp.* have intrinsic elevated MICs to imipenem. MIC results for meropenem, doripenem, and/or ertapenem must be used to determine if these organisms should be forwarded to DCLS for testing.

- Pure suspect CRPA isolates on slant or plate media.
  - 1) Isolate identified as *Pseudomonas aeruginosa*
  - 2) Resistance to at least 1 carbapenem:
    - imipenem, meropenem or doripenem = **MIC of  $\geq 8$  µg/mL**
  - 3) Mucoid *Pseudomonas aeruginosa* isolates will NOT be tested
- Pure suspect CRAB isolates on slant or plate media.
  - 1) Isolate identified as *Acinetobacter baumannii* complex
  - 2) Resistant to at least 1 carbapenem:
    - imipenem, meropenem or doripenem = **MIC of  $\geq 8$  µg/mL**

## • **SPECIMEN TRANSPORT**

- Ship isolates at room temperature.
- Submit a completed DCLS Test Request Form with isolate, to include the following:
  - 1) Complete patient information (including address)
  - 2) Complete submitter information
  - 3) Date of specimen collection
  - 4) Time of specimen collection
  - 5) Specimen source
  - 6) Organism suspected
  - 7) On page 2 of the DCLS test request form check the box for “Carbapenem Resistant Organism” under “Bacteriology ID/ Detection”
  - 8) Outbreak ID number (if applicable)
- Include copy of laboratory’s susceptibility results in specimen package
- **Shipping Address:**  
Division of Consolidated Laboratory Services  
Attn: CRO Team  
600 North 5<sup>th</sup> Street  
Richmond, Virginia 23219

## ● **RESULTS REPORTING**

- Alert value results (**results possibly requiring additional reference laboratory testing**) will be verbally reported within 1 working day of results. An alert notification will be sent to the CDC as well.
- Non-alert-value, positive results (**results NOT requiring additional reference laboratory testing**) will be verbally reported to the submitting laboratory within 2 working days of results.
- Negative results will not be verbally reported.
- A hard copy report of final test results will be provided to the submitter by mail.
  - CRAB isolates will include: identification, AST, and molecular results
  - mCIM positive CRE/CRPA isolates will include: identification, AST, phenotypic, and molecular results
  - mCIM negative CRE/ CRPA isolates will include: identification and phenotypic result only
- VDH will receive electronic reporting of results.
- Reference laboratory results for additional characterization testing will be provided to the submitting laboratory and VDH within 2 working days of receipt from the reference laboratory.

## ● **CONTACT INFORMATION**

- **For questions regarding testing, specimen collection or shipment, please contact one of the following:**
  - Kristin DiBiase (804-648-4480 x: 339)
  - Dr. LaToya Griffin-Thomas (804-648-4480, ext. 281)
  - Microbial Reference Laboratory (804-648-4480, ext. 287 or 251)
  - Molecular Detection and Characterization Laboratory (804-648-4480, ext. 298)
  - DCLS 24/7 Emergency Mobile Number (804-335-4617)

## ● **RESOURCES**

Updated information on CRO testing can be found on the DCLS website, under the “Hot Topics” link:  
<https://dgs.virginia.gov/division-of-consolidated-laboratory-services/updates/hot-topics/>

Link to DCLS Test Request Form:  
<https://dgs.virginia.gov/division-of-consolidated-laboratory-services/resources/submission-forms/>

Link to CDC Antimicrobial Resistance Laboratory Network (ARLN) webpage  
<https://www.cdc.gov/drugresistance/ar-lab-networks/domestic.html>