

Virginia Division of Consolidated Laboratory Services

ANIONS IN WATER BY CHEMICALLY SUPPRESSED ION CHROMATOGRAPHY ASTM D 4327-97					
Facility Name: _____ VELAP ID _____					
Assessor Name: _____ Analyst Name: _____ Inspection Date _____					
Relevant Aspect of Standards	Method Reference	Y	N	N/A	Comments
Records Examined: SOP Number/ Revision/ Date _____ Analyst: _____					
Sample ID: _____ Date of Sample Preparation: _____ Date of Analysis: _____					
Were the nitrite stock solutions refrigerated and prepared monthly?	8.5.5				
Was the sodium sulfate used to make sulfate stocks dried in an oven at 105°C for 1 hour prior to use?	8.5.7				
Were at least three standards and a blank used for calibration?	8.6				
Were working calibration standards prepared fresh daily?	8.6				
When nitrite, nitrate, or phosphate were analytes, were samples refrigerated at 4°C?	9.2				
Was eluent pumped through the system long enough to obtain a stable baseline prior to calibration?	11.2				
Notes/Comments:					

METHOD CHECKLISTS ARE AN INTERVIEW TOOL USED BY ASSESSORS AND ARE NOT TO BE CONSIDERED AS A SUBSTITUTE FOR REQUIREMENTS OF THE PUBLISHED METHOD. CHECKLISTS ARE SUBJECT TO CHANGE.