

| PURGE AND TRAP CAPILLARY-COLUMN GC/MS<br>SM 20 <sup>th</sup> ED 6200  |                           |   |   |     |          |
|---|---------------------------|---|---|-----|----------|
| Facility Name: _____ VELAP ID _____   |                           |   |   |     |          |
| Assessor Name: _____ Analyst Name: _____ Inspection Date _____  |                           |   |   |     |          |
| Relevant Aspect of Standards  | Method Reference          | Y | N | N/A | Comments |
| Records Examined: SOP Number/ Revision/ Date _____ Analyst: _____   |                           |   |   |     |          |
| Sample ID: _____ Date of Sample Preparation: _____ Date of Analysis: _____  |                           |   |   |     |          |
| Were sample containers purchased precleaned or cleaned with detergent, tap, and distilled water and dried for 1-hour at 105°C?                      | 6010 B 1                  |   |   |     |          |
| Were samples collected in duplicate?  | 6010 B 1                  |   |   |     |          |
| Was a replicated sample collected with each sample set?   | 6010 B 1                  |   |   |     |          |
| Were field reagent blanks of reagent water collected and shipped along with sample containers?  | 6010 B 1<br>6200 B 1 b    |   |   |     |          |
| Were samples collected by filling sample containers to just overflowing without trapping air bubbles in sealed bottles?                             | 6010 B 1                  |   |   |     |          |
| Were taps flushed until temperatures stabilized prior to sampling?  | 6010 B 1                  |   |   |     |          |
| Were samples only known to <b>not</b> contain residual chlorine preserved with HCl?   | 6010 B 1 1)               |   |   |     |          |
| Was ascorbic acid used to preserve samples believed to contain residual chlorine where ascorbic acid was known not to cause interferences?          | 6010 B 1 2)               |   |   |     |          |
| Were reagent blanks analyzed to confirm the absences of interferences with each sample batch?   | 6010 B 1 2)<br>6010 B 3 a |   |   |     |          |
| Did reagent blanks include all reagents and preservatives that contact samples, and were they carried through all sample preparations and analyses? | 6010 B 3 a                |   |   |     |          |
| Were samples chilled to 4°C after collection?   | 6010 B 1 2)               |   |   |     |          |
| Were a minimum of 5 concentrations used for calibration?  | 6020 B 1 a<br>6200 B 3 j  |   |   |     |          |
| Notes/Comments:   |                           |   |   |     |          |

| <b>PURGE AND TRAP CAPILLARY-COLUMN GC/MS<br/>SM 20<sup>th</sup> ED 6200</b>                                      |                         |          |          |            |                 |
|--|-------------------------|----------|----------|------------|-----------------|
| <b>Relevant Aspect of Standards</b>  | <b>Method Reference</b> | <b>Y</b> | <b>N</b> | <b>N/A</b> | <b>Comments</b> |
| Did calibration concentrations have no more than one order of magnitude between them?                            | 6020 B 1 a              |          |          |            |                 |
| Was one calibration standard near but above the MDL?   | 6200 B 3 j              |          |          |            |                 |
| Were continuing calibration verification analyzed at least every 20 samples or 12 hours?                         | 6020 B 1 b              |          |          |            |                 |
| Were the acceptance criteria for CCVs between 80 and 120%?   | 6020 B 1 b              |          |          |            |                 |
| Did CCV concentrations vary to confirm the entire calibration range?   | 6020 B 1 b              |          |          |            |                 |
| Were analysis runs closed with an LFB of known concentration?  | 6020 B 1 c              |          |          |            |                 |
| Were LFBs included with each sample batch or set?  | 6020 B 3 b              |          |          |            |                 |
| Were surrogates added before sample preparation so they could be carried through the entire preparation?         | 6020 B 3 d              |          |          |            |                 |
| Were externally prepared QCS analyzed whenever new stock solutions were prepared or every quarter?               | 6020 B 3 e              |          |          |            |                 |
| Were laboratory-fortified samples (LFS) prepared with each sample batch?   | 6020 B 3 f              |          |          |            |                 |
| Were laboratory-fortified sample duplicates prepared from duplicate samples included in each sample batch?       | 6020 B 3 g              |          |          |            |                 |
| Were reagent-blanks analyzed daily, and the values never used to correct sample values?                          | 6200 B 1 b              |          |          |            |                 |
| Were LRBs analyzed after unusually concentrated samples to check for carryover?                                  | 6200 B 1 b              |          |          |            |                 |
| When stock standard compound purities were less than 96%, were the calculated standard concentrations corrected? | 6200 B 3 g              |          |          |            |                 |
| Notes/Comments:  |                         |          |          |            |                 |

| PURGE AND TRAP CAPILLARY-COLUMN GC/MS<br>SM 6200 B 20 <sup>TH</sup>  |                  |   |   |     |          |
|--|------------------|---|---|-----|----------|
| Relevant Aspect of Standards   | Method Reference | Y | N | N/A | Comments |
| Were stock standards stored at -10°C to -20°C with minimal headspace and away from light?  | 6200 B 3 g       |   |   |     |          |
| Were secondary standards prepared weekly for gases?  | 6200 B 3 h       |   |   |     |          |
| Were secondary standards stored with minimal headspace in a freezer?   | 6200 B 3 h       |   |   |     |          |
| Were working calibration standards stored for not longer than 24 hours if in sealed, zero-headspace containers and for not longer than 1 hour if in another type of container? | 6200 B 3 j       |   |   |     |          |
| Was 25 ng BFB analyzed at the beginning of each day before any sample analysis?  | 6200 B 4 b       |   |   |     |          |
| Were all key m/z criteria met in BFB analysis before sample analysis?  | 6200 B 4 b       |   |   |     |          |
| If using internal calibration, were average response factors used only if RSDs were less than 20%?   | 6200 B 4 c 2)    |   |   |     |          |
| If using external calibration, was linearity through zero assumed only if RSDs were less than 20% throughout the entire calibration range?                                     | 6200 B 4 c 3)    |   |   |     |          |
| Were samples brought to room temperature prior to analysis?  | 6200 B 4 d       |   |   |     |          |
| Notes/Comments:  |                  |   |   |     |          |