

Total Residual Chlorine by Iodometric Method II (Iodine Titrant) SM 18th ED 4500-Cl C						Page 1 of 1
Facility Name: _____ VELAP ID: _____						
Assessor Name: _____ Analyst Name: _____ Inspection Date: _____						
Relevant Aspect of Standards	Method Reference	Y	N	N/A	Comments	
<i>Records Examined:</i> SOP Number/ Revision/ Date _____ Analyst: _____						
Sample ID: _____ Date of Sample Preparation: _____ Date of Analysis: _____						
The phenylarsine oxide (PAO) or sodium thiosulfate standardized to 0.00564 N.	4500-Cl C 3.2					
Standard iodine titrate prepared to 0.0282 N and standardized daily prior to use (as described in 4500-Cl C 3.g.).	4500-Cl C 3.h.					
Cl ₂ -demand-free water prepared and stored as specified.	4500-Cl C 3.m.					
Acetate buffer prepared with chlorine-demand-free water.	4500-Cl C 3.e.					
Sample volume for Cl ₂ concentrations < 10 mg/L is 200 mL and for higher concentrations volume is diluted to 200 mL with Cl ₂ -demand-free water so not more than 10 mL of PAO solution is used.	4500-Cl C 4.a.1)					
For titration by iodine titration use 5 mL of PAO or S ₂ O ₃ (Cl ₂ of 2-5 mg/L) or 10 mL (Cl ₂ of 5-10 mg/L)	4500-Cl C 4.a.2)					
Add excess KI (~1 gram)	4500-Cl C 4.a.2)					
Add 4 mL of acetate buffer or enough to reduce the pH to between 3.5 to 4.2.	4500-Cl C 4.a.2)					
Add 1 mL of starch solution.	4500-Cl C 4.b.2)					
Correct equation used for calculation: mg Cl as Cl ₂ /L = ((A-5B) X 200) / C , where: A = mL of 0.00564N PAO or Thiosulfate, B = mL of 0.0282N I ₂ , and C = mL sample.	4500-Cl C.5.a.					
Notes/ Comments:						