

**Dissolved oxygen (D.O.) by Membrane Electrode Method
SM 18th, 19th and 20th 4500-O G**

Page 1 of 1

Facility Name: _____ VELAP ID: _____

Assessor Name: _____ Analyst Name: _____ Inspection Date: _____

Relevant Aspect of Standards**Method
Reference****Y****N****N/A****Comments***Records Examined:* SOP Number/ Revision/ Date _____ Analyst: _____

Sample ID: _____ Date of Sample Preparation: _____ Date of Analysis: _____

If samples are collected from a line under pressure, does the sampler avoid turbulence, formation of bubbles, allow the bottle to overflow two or three times its volume, and stopper it so that no air bubbles are entrained?

4500-O G.1.c,
4500-O B.3

If samples are collected, are samples analyzed immediately (unless known to lack iodine demand)?

4500-O G.1.c,
4500-O B.4

Is the membrane changed frequently?

4500-O G.1.b

Does the analyst take care in changing the membrane to avoid contamination of the sensing element and also trapping of minute air bubbles under the membrane?

4500-O G.3.b

Is the meter calibrated frequently, following manufacturer's calibration procedure exactly?

4500-O G.1.b
4500-O G.3.a

Is there sufficient sample flow across the membrane surface to overcome erratic response?

4500-O G.3.b

Does the analyst follow all precautions recommended by the manufacturer to ensure acceptable results? (This could include factors such as meter stabilization, electrode storage, and type of electrode filling solution used.)

4500-O G.3.b

Notes/ Comments: