

<b>CALCIUM BY SM 18<sup>TH</sup> 3500-Ca D</b> <b>EDTA TITRIMETRIC METHOD</b>					<b>Page 1 of 1</b>
Facility Name: _____ VELAP ID _____					
Assessor Name: _____ Analyst Name: _____ Inspection Date _____					
Relevant Aspect of Standards	Method Reference	Y	N	N/A	Comments
<i>Records Examined:</i> SOP Number/ Revision/ Date _____ Analyst: _____					
Sample ID: _____ Date of Sample Preparation: _____ Date of Analysis: _____					
(NOTE: Polluted waters and wastewater samples may be digested by SM 18 <sup>th</sup> 3030E or SM 18 <sup>th</sup> 3030I prior to analysis by this method.)	3500-Ca D.3.a				
Were hard waters (alkalinity > 300 mg CaCO <sub>3</sub> /L) treated by first diluting OR by neutralizing the alkalinity with acid and then boiling for 1 minute prior to titrating?	3500-Ca D.3.b				
Was sufficient NaOH added to the sample to produce a pH of 12 to 13?	3500-Ca D.3.c				
Was EDTA titrant added slowly with continuous stirring?	3500-Ca D.3.c				
Were samples titrated immediately after adding the alkali (1N NaOH) and Indicator (Either Murexide or Eriochrome Blue Blanck R)?	3500-Ca D.3.b				
When Murexide was used as the indicator, were 1 or 2 excess drops of titrant added to make certain that ne further color change occurred?	3500-Ca D.3.c				
Notes/Comments:					