

**TOTAL PHOSPHORUS BY
LACHAT QUIKCHEM METHOD 10-115-01-1-C
"TOTAL PHOSPHORUS IN KJELDAHL DIGESTS"
COPYRIGHT DATE JULY 13, 1992**

PAGE 1 OF 1

Facility Name: _____ VELAP ID: _____

Assessor Name: _____ Analyst Name: _____ Inspection Date: _____

Relevant Aspect of Standards**Method
Reference****Y****N****N/A****Comments***Records Examined:* SOP Number/ Revision/ Date _____ Analyst: _____

Sample ID: _____ Date of Sample Preparation: _____ Date of Analysis: _____

Were tubes first placed in digester at 160°C for one hour?

Step 4

Was digestion continued for 2.5 hours at 380°C?

Step 5

Was the Ascorbic Acid Reducing Solution prepared fresh every 2 or 3 days?

Reagent 4

Were samples for Nitrate-Nitrite (combined) analysis stored at ≤6°C, preserved to pH<2 with H₂SO₄, and held for no longer than 28 days?CFR 136.3
Table II

Notes/Comments: