

E COLI by Modified m-TEC/ Membrane Filtration HACH 8367 (DOC316.53.001211)					
Facility Name: _____			VELAP ID _____		
Assessor Name: _____		Analyst Name: _____		Inspection Date _____	
Relevant Aspect of Standards	Method Reference	Y	N	N/A	Comments
Records Examined: SOP Number/ Revision/ Date _____ Analyst: _____					
Sample ID: _____ Date of Sample Preparation: _____ Date of Analysis: _____					
1) Are samples refrigerated at <10°C, treated with 0.0008% Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> <sup>5</sup> , and analyzed within 8 hours of collection?	40 CFR 136.3 Table 11				
2) Are membrane filters, sterile, white grid marked, 47mm diameter, with 0.45µm pore size?	Page 4				
3) Are modified m-TEC 47mm prepared agar plates, HACH catalog number 2811815, used?	Page 4				
4) Are sample volumes chosen to produce 20 to 200 colony forming units (CFUs) per filter? ( <i>Ideally, for wastewater, 20 to 80 CFUs.</i> )	Page 1				
5) Is the work bench disinfected with a germicidal cloth, dilute bleach solution, bactericidal spray, or dilute iodine solution, and are hands washed thoroughly with soap and water?	Page 1				
6) Is a membrane filter placed grid side up on the filtration apparatus using sterilized forceps?	Page 2				
7) When the sample is less than 20 mL (diluted or undiluted), is 10mL of sterile dilution water added to the filter funnel before applying vacuum?	Page 1				
8) Prior to filtration, is the sample inverted for 30 seconds to mix?	Page 2				
9) After the sample is filtered, are the funnel walls rinsed three times with 20 to 30mL of sterile buffered dilution water?	Page 2				
10) Is the vacuum released when the filter is dry to prevent damage to the filter?	Page 2				
11) Is the filter aseptically removed and rolled grid side up onto the agar plate to avoid trapping air under the filter?	Page 2				
Notes/ Comments:					

E COLI by Modified m-TEC/ Membrane Filtration HACH 8367 (DOC316.53.001211)					
Facility Name: _____			VELAP ID _____		
Assessor Name: _____		Analyst Name: _____		Inspection Date _____	
Relevant Aspect of Standards	Method Reference	Y	N	N/A	Comments
12) Is the dish inverted and incubated at $35 \pm 0.5^{\circ}\text{C}$ for 2 hours and then at $44.5 \pm 0.2^{\circ}\text{C}$ for 22 hours?	Page 2				
13) After 22 hours, are number of red or magenta colonies counted using a 10 to 15X microscope? ( <i>Red or magenta colonies confirm the presence of E. coli.</i> )	Page 2				
14) Is coliform density calculated using the following equation?  <i>Coliform colonies/ 100 mL =</i> $\frac{\text{colonies counted} \times 100}{\text{mL sample filtered}}$	Page 3				
15) If growth covers the entire filtration area of the membrane, or a portion of it, and colonies are not discrete, are results reported as "Confluent Growth With or Without Coliforms"?	Page 3				
16) If the total number of colonies (coliforms plus non-coliforms) exceeds 200 per membrane or the colonies are too indistinct for accurate counting, are results reported as "Too Numerous to Count" (TNTC)?	Page 3				
17) For either situation above (#15 or #16), is a new sample run using a dilution that will give about 50 coliform colonies and not more than 200 colonies of all types?	Page 3				
18) For nonpotable water, if no filter meets the desired minimum colony count, is the average coliform density calculated?  <i>Coliform colonies/ 100 mL=</i> $\frac{\text{Sum of colonies in all samples} \times 100}{\text{Sum of volumes (mL) of all samples}}$	Page 3				
19) <u>Recommendation:</u> <i>Pseudomonas aeruginosa</i> is recommended as a negative control and <i>Escherichia coli</i> as a positive control.	Page 3				
Notes/Comments:					