

NITRATE ELECTRODE METHOD SM 18th and 19th 4500-NO₃⁻ D						Page 1 of 1
Facility Name: _____ VELAP ID _____						
Assessor Name: _____ Analyst Name: _____ Inspection Date _____						
Relevant Aspect of Standards	Method Reference	Y	N	N/A	Comments	
<i>Records Examined:</i> SOP Number/ Revision/ Date _____ Analyst: _____						
Sample ID: _____ Date of Sample Preparation: _____ Date of Analysis: _____						
Were unpreserved samples stored as follows? <i>Nonpotable: ≤ 6°C up to 48 hours</i> <i>Drinking water: 4°C up to 48 hours unless chlorinated, which can be held up to 14 days</i>	40CFR 136.3, 40CFR 141.23					
Were samples held longer than 24 hours preserved with H ₂ SO ₄ , stored at ≤6°C, and analyzed within 28 days?	40CFR 136.3, 40CFR 141.23					
When NO ₃ ⁻ and NO ₂ ⁻ were determined as separate species, were samples never acidified?	4500-NO ₃ A 2 Introduction					
Did buffer solution contain Ag ₂ SO ₄ , sulfamic acid, Al ₂ (SO ₄) ₃ , and a buffer at pH 3?	4500-NO ₃ D 1 b					
Was the pH meter capable of 0.1 mV resolution?	4500-NO ₃ D 2 a					
Was the outer chamber of the reference electrode filled with (NH ₄) ₂ SO ₄ ?	4500-NO ₃ D 2 b					
Were the manufacturer's instructions followed for the care and storage of the nitrate ion electrode?	4500-NO ₃ D 2 c					
Were electrode tips allowed to stabilize in samples and solutions about 1 minute prior to the recording of millivolt readings?	4500-NO ₃ D 4 a/b					
Were electrode tips rinsed and blotted dry between samples?	4500-NO ₃ D 4 a					
Were calibration curves made by semilogarithmic plots with concentration on the logarithmic axis and millivolt responses on the linear axis?	4500-NO ₃ D 4 a					
Were calibration slopes +57 ± 3 mV/decade?	4500-NO ₃ D 4 a					
Notes/Comments:						